

Product Name :
Instrument Trolley**Product Code :**
NMTOY10003**Description :**

Instrument Trolley Frame work made of ERW tubular pipe mounted on four 7.5 cms. dia. castors.
Two fixed stainless steel shelves of heavy gauge.
Pre-treated and epoxy powder coated.

Also available :

Instrument Trolley - All stainless steel (NMTOY10002A-0086)

Instrument Trolley Without railing. (NMTOY10002B-0087)

Instrument Trolley Size: 68L x 45W x 82H cms.

Instrument Trolley Size: 76L x 45W x 82H cms.

Instrument Trolley Size: 91L x 45W x 82H cms.