

Product Name :
Dressing Trolley**Product Code :**
NMTOY10002**Description :**

Dressing Trolley: Tubular pipe frame work. Two shelves if stainless steel with railings. Fitted with S.S. basin and S.S. pail. Mounted on castor wheels.

Dressing Trolley Steel: Stainless steel tubular frame work mounted on four 100 mm dia swiveling castors. Two stainless steel shelves with four side railings. Supplied with stainless steel bowl and bucket.

NAUGRA[®]**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,