

NAUGRA[®]

Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Serum Test Tubes (MOQ: 100,000 Pcs Assorted Size)

Product Code :
NMCBS35084



Description :

Serum Test Tubes (MOQ: 100,000 Pcs Assorted Size) NET Plain tubes with clot activators are used for collecting blood samples and extracting serum by centrifugation, for clinical examinations

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: [SITE](#)