



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
ESR Vacuum Blood Tubes (MOQ: 100,000 Pcs Assorted Size)

Product Code :
NMCBS35089



Description :

ESR Vacuum Blood Tubes (MOQ: 100,000 Pcs Assorted Size) ESR Vacuum Blood Tubes (MOQ: 100,000 Pcs Assorted Size) ESR is a simple, non-specific screening test utilized for monitoring the presence of inflammation in the body. ESR tubes contain tri-sodium citrate solution and are used for blood sedimentation analysis. Mixing recommendation: ESR tube should be gently inverted 180° back and forth 8-10 times



Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,