

NAUGRA[®]

Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Safety Box - Cardboard, 5 Litre

Product Code :
NMCBS35105



Description :

Safety Box - Cardboard, 5 Litre Volume : 5.0Litres safety box Capacity : 150 auto disable syringes of 0.5ml with needles Wall thickness : 1.4mm, Size : 285x160x130mm(+/-2mm) Printing : 2 colors standard printing. white background, printing red and black Diameter of syringe insert hole : 35-40mm Material : Recycled board, water resistant with waterproof membrane

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,