

NAUGRA[®]

Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Reusable Reinforced Silicone Laryngeal Mask Airway

Product Code :
NMCBS35135



Description :

Reusable Reinforced Silicone Laryngeal Mask Airway 100% Medical Grade Silicone. A Stainless Steel spring is inserted in the tube to reduce the risk of kinking. Can be autoclaved repeatedly at 134°C. Reusable up to 40 times

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: [SITE](#)