

Product Name :
Semi Fowler Steel Bed Manual**Product Code :**
NMHOPB00029**Description :**

Semi Fowler Steel Bed Manual India China

Technical Specification :

Semi Fowler Steel Bed Manual with mechanically operated back rest by super-smooth crank mechanism. Size: L 2120mm

W 940mm

H 500mm. Removable & interchangeable high quality beautiful ABS engineering plastic head panel and foot panel. Epoxy coated mild steel frame work and 2 section top. Mounted on protective stumps. Provision for I.V. Rod on both sides of the bed. Freight saving knock-down construction