

**Product Name :**  
Semi Fowler Steel Bed Manual**Product Code :**  
NMHOPB00029**Description :**

Semi Fowler Steel Bed Manual India China

**Technical Specification :**

Semi Fowler Steel Bed Manual with mechanically operated back rest by super-smooth crank mechanism. Size: L 2120mm

W 940mm

H 500mm. Removable &amp; interchangeable high quality beautiful ABS engineering plastic head panel and foot panel. Epoxy coated mild steel frame work and 2 section top. Mounted on protective stumps. Provision for I.V. Rod on both sides of the bed. Freight saving knock-down construction