

**Product Name :**  
Standard Semi-Fowler Bed**Product Code :**  
NMHOPB00031**Description :**

Standard Semi-Fowler Bed Standard semi-fowler bed with mechanically operated back rest by crank mechanism. Size: 2030 ĩ½ 900 ĩ½ 530 mm. Epoxy coated mild steel frame work and 2 section top. Tubular head & foot bows made of stainless steel. Mounted on protective stumps. Provision for I.V. Rod. Freight saving knock-down construction

**NAUGRA**<sup>®</sup>**Medical Hospital Equipments Manufacturers**Website: [SITE](#)