

**Product Name :**  
Naugramedical steel Fowler Bed Manual**Product Code :**  
NMHOPB00042**Description :****Naugramedical steel Fowler Bed Manual India China****Technical Specification :**

Naugramedical steel Fowler Bed Manual with mechanically operated back rest and knee rest by super-smooth crank mechanism. Size: L 2140mm W 940mm

H 530mm. Four easy lifting guardrails (2 on each side), which are safe and reliable, and can be fixed upward and downward. Removable & interchangeable high quality beautiful ABS engineering plastic head panel and foot panel. Epoxy coated mild steel frame work and 4 section perforated top. 125mm dia noiseless castors with Central braking system. Provision for I.V. Rod