

Product Name :
Fowler Provision Bed Manual 2 Function

Product Code :
NMHOPB00045



Description :

Fowler Provision Bed Manual 2 Function

Technical Specification :

Fowler Provision Bed Manual 2 Function with mechanically operated back rest and knee rest by super-smooth crank mechanism. Size: L 2140

W 940

H 500mm. Removable & interchangeable high quality beautiful flat tube head panel and foot panel. Collapsible side railings. (Top bar Aluminium, Vertical bars SS, Bottom pipe epoxy) Epoxy coated mild steel frame work and 4 section perforated top. 125mm dia noiseless castors with Central braking system. Provision for I.V. Rod on both sides of the bed