

**Product Name :**  
Fowler Bed Manual Steel Two Function**Product Code :**  
NMHOPB00049**Description :**

Fowler Bed Manual Steel Two Function

**Technical Specification :**

Fowler Bed Manual Steel Two Function with mechanically operated back rest and knee rest by super-smooth crank mechanism. Size: L 2120mm

W 980mm

H 500mm. Removable &amp; interchangeable high quality beautiful ABS engineering plastic head panel and foot panel. Collapsible side railings. (Top bar Aluminium, Vertical bars SS, Bottom pipe epoxy) Epoxy coated mild steel frame work and 4 section top. 125mm dia noiseless castors (all with individual brakes) Provision for I.V. Rod on both sides of the bed. Freight saving knock-down construction