

**Product Name :**  
Fowler Bed, Manual (with Provision for Bed Pan)**Product Code :**  
NMHOPB00062**Description :**

Fowler Bed, Manual (with Provision for Bed Pan) Fowler Bed with mechanically operated back rest 0 to 75°½ and knee rest 0 to 40°½ by crank mechanism. Size: 2000 x 900 x 500 mm. Provision for Bed pan. With side railings. Epoxy coated mild steel frame work and 4 section top. Mounted on protective stumps. Provision for I.V. Rod on both sides of the bed. Freight saving knock-down construction