

Product Name :
Orthopaedic Bed, Fowler**Product Code :**
NMHOPB00068**Description :**

Orthopaedic Bed, Fowler Orthopaedic bed with mechanically operated back rest & knee rest all by Super Smooth crank mechanism. Size: 2150 i_l½ 940 i_l½ 500 mm. Epoxy Coated Steel frame work and 4 section perforated top. Stainless Steel Orthopaedic attachment (without weight chain and weight hangers). Head panel and foot panel equipped with safety lock and roller bumpers. Mounted on protective stumps. Provision for I.V. Rod. Freight saving knock-down construction