

NAUGRA[®]

Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Hospital Doctor Room Chair

Product Code :
NMCBS10003



Description :

Hospital Doctor Room Chair M.S. tubular frame with heavy polypropylene base mounted on 5 cms dia. castors. With Cushioned seat and back. Lever controlled independant seat tilt and back angle. Pneumatic height adjustable. Pre-treated and epoxy powder coated. Also Available in black leather or a number of fabric choices. Hospital Doctor Room Chair, Hospital Chairs and Stools Exporters, Hospital Chairs and Stools Manufacturers

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: [SITE](#)