

**NAUGRA**<sup>®</sup>

Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Infant Bed / Child Cot / Baby Bassinet

**Product Code :**  
BCIVP10002



**Description :**

Infant Bed / Child Cot / Baby Bassinet Size: L 105mm  $\pm$  1/2 W 45mm  $\pm$  1/2 H 1045mm  $\pm$  1/2 855mm (Height Adjustable by hydraulic system). In build touch keyboard. Weight scale to Measure baby weight . Tilt angle 15 $\pm$  1/2, easy for mother and baby to communicate and feed. Hand rail could be placed flat or at 90 $\pm$  1/2 angle. Can be used as a mosquito net support. Additional Features: SD memory card to play music, can store music in MP3 format Atmosphere temperature measuring

**NAUGRA**<sup>®</sup>

**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,