

NAUGRA[®]

Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Infant Bed / Child Cot / Baby Bassinet

Product Code :
BCIVP10002



Description :

Infant Bed / Child Cot / Baby Bassinet Size: L 105mm \pm 1/2 W 45mm \pm 1/2 H 1045mm \pm 1/2 855mm (Height Adjustable by hydraulic system). In build touch keyboard. Weight scale to Measure baby weight . Tilt angle 15 \pm 1/2, easy for mother and baby to communicate and feed. Hand rail could be placed flat or at 90 \pm 1/2 angle. Can be used as a mosquito net support. Additional Features: SD memory card to play music, can store music in MP3 format Atmosphere temperature measuring

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: [SITE](#)