

NAUGRA[®]

Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Infant Bed / Child Cot S.S Adjustable Heights

Product Code :
BCIVP10006



Description :

Infant Bed / Child Cot S.S Adjustable Height Size: 750 ĩ½ 400 ĩ½ 850/1000 mm. Height adjustable approx. 850 - 1000 mm. Completely made of stainless steel. 4 swivel castor (all with brakes)

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: [SITE,](#)