

Product Name :
Semi-Fowler Bed For Children, with Side Railings**Product Code :**
BCIVP10008**Description :**

Semi-Fowler Bed For Children, with Side Railings Semi-Fowler Bed with mechanically operated Back rest by crank mechanism. Size: 1880 x 1/2 880 x 1/2 580 mm. Epoxy coated mild steel frame work with 2 section top. S.S. side railings with Up & Down movement. With S.S. Shelf at the bottom for keeping blankets, basin etc. Tubular Head & Foot Bows also made of stainless steel. Mounted on four silent swivel castor (all with brakes). Provision for I.V. Rod. Freight saving knock-down construction

NAUGRA[®]**Medical Hospital Equipments Manufacturers**Website: SITE