

NAUGRA[®]

Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Overbed Table (Adjustable With Geared Handle)

Product Code :
NMIVP20031



Description :

Overbed Table (Adjustable With Geared Handle) Table base and frame work made of strong rectangular / square steel tubes mounted on four 5 cms dia. castors. Overall size of laminated top : 76L x 40W cms. adjustable by gear handle. Pre-treated and epoxy powder coated. Knockdown construction

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: [SITE](#)