

NAUGRA[®]

Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Overbed table, Height Adjustable

Product Code :
NMIVP20033



Description :

Overbed table, Height Adjustable Overbed table with Laminated Wooden Top. Size: 795 ÷ 1/2 400 ÷ 1/2 710 / 970 mm. Framework is made of Epoxy steel. Adjustable height from 710 - 970 mm by gas spring. 4 Swivel castor

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: [SITE](#)