

Product Name :
Ambulance Foldaway Stretcher DX**Product Code :**
NMSTR10003**Description :**

Ambulance Foldaway Stretcher DX with wheels, outriggers and can be folded. It is made of high strength aluminum alloy and the water-tightness leather material. The back of this stretcher can be lifted, and the maximum angle is 60° ± 1/2°. It is characterized by its being fashionable structure, light-weighted, easy-carried, use-safely, reliable quality and easy for sterilization.

Dimensions

Product size (L x W x H): 185 x 50 x 24cm

Folding size (L x W x H): 92.5 x 50 x 11cm

Packing size (1pcs): 96 x 54 x 12cm Ambulance Foldaway Stretcher DX, Hospital Stretchers Exporters, Hospital Stretchers Manufacturers