

Product Name :
Folding Stretcher**Product Code :**
NMSTR10006**Description :**

Folding Stretcher with outrigger is made of high strength aluminum alloy material and the waterproof leather materials. This stretcher with outrigger can be folded in length and width, easy for transportation and storage. The orange fabric is fire proof, water proof, anti-cracking and easy to clean. It is characterized by its being light weight, portable, reliable quality, easy to carry out and safe to use, etc. It is mainly used for the hospitals, emergency centers and battle field to carry patients and wounded person.

Dimensions:

Product size (L x W x H): 208 x 54.5 x 13cm

Folding size (L x W x H): 104 x 9 x 17cm

Packing size (2pcs): 108 x 18 x 18.5cm

N.W.: 4.7kg

G.W.: 5.3kg

Load bearing: less than 159kg Folding Stretcher, Hospital Stretchers Exporters, Hospital Stretchers Manufacturers

