

**Product Name :**  
Manual Bed, Two Function For Home Care**Product Code :**  
NMHOPB00075**Description :**

Manual Bed, Two Function For Home Care Two function manual bed for home care with crank operated back rest

**Technical Specification :**

- Two function manual bed for home care with crank operated back rest:0°-70°±2°, knee rest:0°-40°±2°
- Size: 2080mm(L)\*1012mm(W)\*500mm(H)
- Epoxy coated mild steel frame work and 4 section perforated top.
- Collapsible side railings. (Top bar Aluminium, Vertical bars SS, Bottom pipe epoxy)
- Removable high quality beautiful wooden head panel and foot panel.
- 125mm dia noiseless castor all with brakes.
- Provision for I.V. Rod on all sides of the bed.