



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Professional Ophthalmoscope

Product Code :
NMDIG10047



Description :

Professional Ophthalmoscope Halogen light for true tissue color and consistent, long-lasting illumination.
Apertures: large and small spot, semicircle, and red-free filter, Correcting lenses range: +in 1D steps
|1-4|6|8|10|15|20|0 & -in 1D steps |1-4|6|8|10|15|20 Bulb:2.5V 1.6W Halogen Lamp



Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,