

NAUGRA[®]

Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Fingertip Pulse Oximeter

Product Code :
NMDIG10049



Description :

Fingertip Pulse Oximeter Dependable Oxygen Level and Pulse Rate measurements Small & lightweight design Low power consumption Large LED Display of SpO₂, PR, and Pulse Bar Auto shut-off with no signal for 5 seconds Suitable for adult, pediatric and neonate

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,