

**NAUGRA**<sup>®</sup>

Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Fingertip Pulse Oximeter

**Product Code :**  
NMDIG10050



**Description :**

Fingertip Pulse Oximeter Dependable Oxygen Level and Pulse Rate measurements LCD screen display SpO2, PR, and SpO2 waveform Small & lightweight design Low power consumption Auto shut-off with no signal for 5 seconds Build-in memory, with PC software analyzing Suitable for adult, pediatric and neonate

**NAUGRA**<sup>®</sup>

**Medical Hospital Equipments Manufacturers**

Website: [SITE](#)