

Product Name :
Fetal Doppler**Product Code :**
NMDIG10053**Description :**

Fetal Doppler Doppler can meet routine examination requirement of obstetricians, midwives, and pregnant women. It is a high performance model with FHR digital LCD display. It has three heart rate processing modes: real-time FHR display mode, averaged FHR display mode and manual mode. Both of the models have audio output, and can connect with earphone or recorder with audio input. It has a slim, power saving and interchangeable ultrasound probe

Feature: Easy to use Accurate FHR detection and LCD display Auto and manual counting Built-in speaker Automatic switch-off Battery energy indication Function: FHR Detection: Speaker output/Earphone output LCD display, Three working modes: Real-time FHR Display mode, Averaged FHR Display mode, Manual mode, continuous operation time?15hours

NAUGRA[®]**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,