

NAUGRA[®]

Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Fetal Doppler

Product Code :
NMDIG10054



Description :

Fetal Doppler Doppler can meet routine examination requirement of obstetricians, midwives, and pregnant women. It has manual mode heart rate processing mode. It has audio output, and can connect with earphone or recorder with audio input. It has a slim, power saving and interchangeable ultrasound probe. Feature: Easy to use Auto and manual counting Built-in speaker Automatic switch-off Battery energy indication Function: Detection: Speaker output/Earphone output, Manual mode, continuous operation time?15hours

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,