

**Product Name :**  
12 Lead Handheld Single Channel ECG Machine**Product Code :**  
NMDIG10057**Description :**

12 Lead Handheld Single Channel ECG Machine

**Technical Specification :****12 Lead Handheld Single Channel ECG Machine**

Input circuit: Floating  
Calibrating Voltage :1mV  
Input current: 103dB  
A/D: 12bit  
Frequency Response: 0.05Hz–150Hz  
Filter: AC 50/60Hz EMG: 35-45Hz  
Skin Voltage Tolerance: +300mV  
Sensitivity Setting 5, 10, 20mm/mV, conversion deviation 25%  
Noisy Level: ?15μVp-p  
Safety standard: IEC class I, type CF  
Lead selection: Automatic and Manual  
Operation mode: Automatic and Manual  
Display: LCD 160\*96 visual size: 34\*60  
Recording system: Thermal-array printing  
Power supply: AC220V±10%  
DC 7.4 lithium battery

---

**NAUGRA<sup>®</sup>**

**Medical Hospital Equipments Manufacturers**

Website: [SITE](#)