

Product Name :
12 Lead Handheld Single Channel ECG Machine

Product Code :
NMDIG10057



Description :

12 Lead Handheld Single Channel ECG Machine

Technical Specification :

12 Lead Handheld Single Channel ECG Machine

Input circuit: Floating
Calibrating Voltage :1mV
Input current: 103dB
A/D: 12bit
Frequency Response: 0.05Hz–150Hz
Filter: AC 50/60Hz EMG: 35-45Hz
Skin Voltage Tolerance: +300mV
Sensitivity Setting 5, 10, 20mm/mV, conversion deviation 25%
Noisy Level: ?15μVp-p
Safety standard: IEC class I, type CF
Lead selection: Automatic and Manual
Operation mode: Automatic and Manual
Display: LCD 160*96 visual size: 34*60
Recording system: Thermal-array printing
Power supply: AC220V±10%
DC 7.4 lithium battery

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,