

**Product Name :**  
Stretcher on Trolley with Side Railings**Product Code :**  
NMSTR10008**Description :**

Welded tubular steel framework. Dished sheet steel top with propelling handles at each end. Mounted on 15cm dia. Castor wheels Upholstered. Supplied with 35mm thick PU mattress. With side railings.

Available in : EPC Finished with S.S. Top, All EPC, All Stainless Steel Stretcher on Trolley with Side Railings, Hospital Stretchers Exporters, Hospital Stretchers Manufacturers

**NAUGRA**<sup>®</sup>**Medical Hospital Equipments Manufacturers**Website: [SITE](#)