

Product Name :
Ultrasonic Nebulizers

Product Code :
NMDIG10101



Description :

Ultrasonic Nebulizers Durable, stretchable double tubes can produce a dense and tiny-particle mist. Power Voltage : AC220V \pm 10% Ultrasonic frequency : 1.7MHz \pm 10% Power : ? 60W Maximum nebulizing rate : ? 3ml/min. Timing range : 0-60 min. no level for adjustment. Continuous loading time : ? 4h Max medication cup capacity : 150ml

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,