

NAUGRA[®]

Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Ultrasonic Nebulizers

Product Code :
NMDIG10102



Description :

Ultrasonic Nebulizers Durable, stretchable tube can produce a dense and tiny-particle mist. Power Voltage : AC220V \pm 10% Ultrasonic frequency : 1.7MHz \pm 10% Continuous working time : ? 4h Consumption power : ? 50W Max medication capacity for small cup : 150ml Max medication capacity for big cup : 350ml Max nebulizing rate: ? 3ml/min

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,