

NAUGRA[®]

Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Air Compressor Nebulizer

Product Code :
NMDIG10103



Description :

Air Compressor Nebulizer This equipment sprays liquid medication to fog panel by compressing air, and forms tiny particles, which flow into the throat through imbibing tube. It adopts low noise, oil-free-high-efficiency valve with no need for daily lubrication forms. It is an ideal product for family and medical unit. Use Technical Data
Maximum compressor pressure : ? 0.15Mpa Compressor free air flow : ? 10L/min Maximum nebulization rate : ? 0.20ml/min Sound level : ? 65db(A)

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: [SITE](#)