

NAUGRA[®]

Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Digital Clinical Thermometer

Product Code :
NMDIG10112



Description :

Digital Clinical Thermometer Numeric display - easy to read With Beep Alert Mercury Free (environmental friendly) Unbreakable High Precision Measurement Range: 32-42½°C & 89.6-109.4½°F Display Resolution: 0.1½°C & 0.1½°F Easy Replaceable Battery Auto power off 8 minutes Last measurement memory In plastic case Digital thermometer can be used orally, rectally & under armpit

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,