

NAUGRA[®]

Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Forehead Thermometer

Product Code :
NMDIG10113



Description :

Forehead Thermometer Display range: 34.0~44.0oC/93.2~111.2oC Accuracy: $\pm 0.3oC/\pm 0.5oF$ Min. scale: 0.1 Memory: 10 memories Measurement time (reference only): 1 second Beeper function Auto shut-off (1min after idle) Battery: 3V, Button battery(LR/CR 2032) Battery life: approx. 4,000 takes Dimensions: 140 x35x 35mm N.W.: 70g (including batteries) Safe and hygienic

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,