

NAUGRA[®]

Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Surgical Operation Light 20

Product Code :
NMOTL10004



Description :

Dome Diameter : 20in
Light intensity : 90000 Lux
Halogen Bulb : 24V - 150W
Reflector : Diachronic
Colour Temp. : 4300K
Temp. Increased at head area : 2~4C
Finish : Oven Baked Surgical Operation Light 20, Hospital OT Lights & Lamps Exporters, Hospital OT Lights & Lamps Manufacturers

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: [SITE](#),