

**NAUGRA**<sup>®</sup>

Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Mobile OT Light

**Product Code :**  
NMOTL10005



**Description :**

Intensity: 70,000 Lux  
Halogen Bulb: 24V - 70W  
Reflector: Diachronic  
Colour Temp.: 4300K  
Temp. Increase at Head Area: 2-4C  
Doma Dia: 22in  
Finish: Oven Baked Mobile OT Light, Hospital OT Lights & Lamps Exporters, Hospital OT Lights & Lamps Manufacturers

**NAUGRA**<sup>®</sup>

**Medical Hospital Equipments Manufacturers**

Website: [SITE](#),