

**Product Name :**  
Ophthalmoscope with Otoscope**Product Code :**  
NMDIG10137**Description :**

Ophthalmoscope with Otoscope Halogen Deluxe Ophthalmoscope 2.5V Xenon Halogen Light for true tissue color and consistent, long lasting illumination 24 color coded Lenses 25 to + 40 for excellent resolution (Red for - Green for +) Rubber eye brow rest to prevent scratching of eye-glasses illumination lens dial for clear identification of adopter setting. Supply with Otoscope handle for 2 x C Size Battery, 2.5v vacuum light (German Bulbs) High quality magnifying lens and 3 specula (size 2.5, 3.5 and 4.5) Supply : Black soft pouch or hard plastic case

**NAUGRA**<sup>®</sup>**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,