

Product Name : Artificial Resuscitator (Ambu Type Bag), Silicone, Autoclavable - Deluxe Quality (Adult)	Product Code : NMANI10005
----------------------------------------------------------------------------------------------------------------------	-------------------------------------



Description :
Artificial Resuscitator (Ambu Type Bag), Silicone, Autoclavable - Deluxe Quality (Adult) Ambu Type Bag Capacity 1600ml. With translucent face mask size: 4. 100% latex free
Technical Specification :
Resuscitator Bags, which is easier to hold, easier to use, completely autoclave for reuse.
These are available in three sizes such as for manual ventilation of adults, children and infants.