

<b>Product Name :</b> Artificial Resuscitator (Ambu Type Bag), Silicone, Autoclavable - Deluxe Quality (Adult)	<b>Product Code :</b> NMANI10005
--	-------------------------------------



<b>Description :</b>
Artificial Resuscitator (Ambu Type Bag), Silicone, Autoclavable - Deluxe Quality (Adult) Ambu Type Bag Capacity 1600ml. With translucent face mask size: 4. 100% latex free
<b>Technical Specification :</b>
Resuscitator Bags, which is easier to hold, easier to use, completely autoclave for reuse.
These are available in three sizes such as for manual ventilation of adults, children and infants.