

Product Name :
Face Masks With S.S. Ring

Product Code :
NMANI10013



Description :

Face Masks With S.S. Ring Size 1. Black Rubber. Padded. CE certified

Technical Specification :

Anesthesia mask. Medical mouth piece mask. Rubber. 5 point metal ring for easy and tight fixation with head straps.

Ideal for use with all types of automatic and manual resuscitators or for anaesthesia.

Ergonomically modelled for a perfect fit on the patient's face and optimized grip for the operator.

With inflatable cushion to ensure perfect fit around the mouth.