



Email : sales@naugramedical.com

Phone: 9896600003

Product Name :
Ceiling Surgical Operating Light Small

Product Code :
NMOTL10008



Description :

Ceiling Surgical Operating Light Small Operation Reflector Lamp offers the illumination as compared to other standard models. Halogen Lights suitable for Major as well as Minor Surgery, fully spring balanced cool lights. Sterilizable handle has been provided in the middle of the dome for adjust the focus. A Shadow less light is suitable for all sizes of Operation Theatres as it gives maximum light output with minimum heat dissipation. Parabolic reflectors & color correcting filters.

Dome Head: 550 mm x 2

Halogen Bulbs: 24V x 70 W x 8 nos.

Light intensity at 1 mt: 1,00,000 Lux

Height Adjustment: 914 mm

Parabolic Reflector: 180 mm

Color Temperature: 4200 K

Control Unit: Transformer Ceiling Surgical Operating Light Small, Hospital OT Lights & Lamps Exporters, Hospital OT Lights & Lamps Manufacturers



Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,

