

**Product Name :**  
Face Masks, Autoclavable, Silicon For Deluxe Quality

**Product Code :**  
NMANI10028



**Description :**

Face Masks, Autoclavable, Silicon For Deluxe Quality Deluxe. Size 4. CE certified

**Technical Specification :**

Lightweight as well as ideally designed and shaped face masks which feature transparent hood that allows easy monitoring of foreign body, vomit, blood excreted through mouth region as well as allow anesthetist to monitor the color of the lips stating the health of the patient.

Further, the strong clear polycarbonate hood also provides for comfortable usage.