

**Product Name :**  
Reservoir Bag, PVC with Adapter

**Product Code :**  
NMANI10031



**Description :**

Reservoir Bag, PVC with Adapter To fit our ambu type bags codes: AN002c, AN002f, AN002g & AN002h

**Technical Specification :**

Non-Rebreather Mask with Reservoir-Adult Size  
Low-resistance check valves prevent rebreathing and allow exhaled gas to escape  
Safety vent allows entrainment of room air  
Clear, soft vinyl for patient comfort and visual patient assessment