

**Product Name :**  
Rebreathing Bags - Black Rubber

**Product Code :**  
NMANI10034



**Description :**

Rebreathing Bags - Black Rubber

**Technical Specification :**

Breathing bag(rebrathing bag) comes in Black color and 100% latex free.

During the process we use high quality medical grade rubber.

These ranges of rebreathing Bag are highly demanded by all medical user, hospital and supplier.

These Anesthesia Rebreathing Bags, act like a reservoir for expired gases and are used to push gases in the lungs.

It help to Reducing the risk of trauma, Easy cleaning and safe usages.