

**Product Name :**  
Ambu Type Bag (Artificial Resuscitator), Adult  
(Autoclavable) Black Rubber

**Product Code :**  
NMANI10036



**Description :**

Ambu Type Bag (Artificial Resuscitator), Adult (Autoclavable) Black Rubber Autoclavable. Optional : Reservoir bag - AN403 (Order separately). CE certified

**Technical Specification :**

Rubber Ambu Bag (Resuscitator) kit consists of manual Resuscitator, Tubing, Reservoir Bag & Mask.

It provides excellent bag re-expansion property. Pressure limiting valve ensures crucial adjustments for enhanced performance.

It is available in two sizes : Adult & Child.