

Product Name : Ambu Type Bag (Artificial Resuscitator), Adult, Black Rubber	Product Code : NMANI10038
---	-------------------------------------



Description :
Ambu Type Bag (Artificial Resuscitator), Adult, Black Rubber Black rubber. Optional : Reservoir bag - AN403 (Order separately). CE certified
Technical Specification :
Rubber Ambu Bag (Resuscitator) kit consists of manual Resuscitator, Tubing, Reservoir Bag & Mask. It provides excellent bag re-expansion property. Pressure limiting valve ensures crucial adjustments for enhanced performance. It is available in two sizes : Adult & Child.