

**Product Name :**  
Ceiling OT Reflector Single**Product Code :**  
NMOTL10010**Description :**

Ceiling OT Reflector Single provides the superior performance and is good for all kinds of Surgery. Inbuilt is an Imported Dichoric Glass Reflector, which increases the Lux up to 80,000. Lights arms are made of Stainless Steel Fitting. During the long time operation, there is no fatigued feeling to the eyes, due to the range of the color temp - 4,200 K as natural color tone Sterilizable handle has been provided in the middle of the dome for regulating the diameter of focus. Light intensity can be adjusted by the Digital Control Panel. Diffuser spreads the light uniformly. It protects the light and gives excellent shadow reduction. It also reduces infrared radiation and enhances the light.

Dome Head: 515 mm

Halogen Bulbs: 24V x 70 W x 3 nos

Light intensity at 1 mt: 80,000 Lux

Height Adjustment: 1014 mm

Dichoric Glass Reflector: 200 mm

Color Temperature: 4200 K

Control Unit: SMPS Ceiling OT Reflector Single, Hospital OT Lights &amp; Lamps Exporters, Hospital OT Lights &amp; Lamps Manufacturers

---

Website: [SITE](#)

**NAUGRA<sup>®</sup>**

**Medical Hospital Equipments Manufacturers**